



### PROSTHETIC KIT

Complete line with drivers, connections and torque wrench.



20.503

### SURGICAL GUIDES KIT

Assists in diagnosing biological spaces between teeth and implants.



20.504

### BONE GRAFT KIT OR ANCLAJE ORTODÓNICO KIT

For bone grafting surgeries and/or orthodontic anchorage.



20.505

20.507

For 22 years, Signo Vinces® has invested in technology and the development of innovative products, to provide its customers with a new way of doing their work, saving time, guaranteeing results and a vision of the future. With you, Signo is much more!



### CERTIFICATIONS

RDC 665/2022 | ANVISA

Boas Práticas de Fabricação de Produtos para a Saúde

EN ISO 13485:2016

Medical Devices - Quality Management Systems Requirements for Regulatory Purposes



# BE DIGITAL



20.200

### BONE EXPANDER KIT

Autoclavable case; Basic surgical kit for different split crest expansions and atraumatic sinus lift.

### IMPLANT REMOVAL KIT

Autoclavable anodized aluminium case. Complete toolkit for removal of fractured implants and screws.

20.320



MATRIZ: SIGNO VINCES EQUIPAMENTOS ODONTOLÓGICOS LTDA  
0800 600 5999 | +55 41 3032.5999  
Rod. Raul Azevedo de Macedo, 733 | Campo Largo / PR | Brasil | 83606-482  
signovinces@signovinces.com.br | www.signovinces.com.br

FILIAL: SIGNO VINCES EUROPA, LDA  
LISBOA | PORTUGAL: Tel +351 217 975 655 | Cel +351 935 014 972  
Av. 5 de Outubro, 151 | 6B | 1050.053  
MADRID | ESPANHA: Teléfono +34 910 845 275 | Teléfono Móvil +34 610 199 040  
Calle Zeus, 6, 2º piso - Oficina 207  
Codigo Postal 28341 Valdemoro | Madrid | Espanha  
europa@signovinces.com.br

84.205 | Rev. 00/2023

DIGITAL LIBRARY: WWW.SIGNOVINCES.COM.BR

# KITS EVO®

## DESIGN, INNOVATION AND PRACTICALITY

With a unique design and modern materials, the kits provide lightness and agility in all procedures. Practical, compact and functional, intuitive.

### EVO GUIDE KIT

For guided surgery of Signo Vinces® implants: Infra, Ares, Duo and Duocon, of measures 3.8 and 4.6mm, lengths 8.5 | 10.0 | 11.5 | 13.0mm.

A simple and compact kit, used for less invasive surgeries than conventional ones, where the virtual planning through radiographic images, allows precise incisions, with the help of a guide previously planned and printed 3D guide.



20.509

NEW


### SURGICAL KIT

Drills and instruments for all Signo Vinces® implants.


A rational and versatile kit, the conical drill set can be used for all conical design implants and the cylindrical drill set can be used in all cylindrical design implants.



20.502




### DIRECT SEQUENCE



### INTERMEDIATE SEQUENCE



### DIRECT SEQUENCE



### INTERMEDIATE SEQUENCE

